

Practical 7 ST

anaerobic cultivation

- Bacteria in relation to oxygen
- Anaerobic infections - sampling and transport
- Cultivation
 - establishing of anaerobic environment
 - anaerobic culture media,
 - procedure: inoculation, isolation

Anaerobes

- are the predominant components of the skin's and mucous membranes normal flora
 - cause infections of endogenous origin
- mixed infections caused by numerous aerobic and anaerobic bacteria are often observed in clinical situations.

Anaerobic infections

- Usually endogenous bacteria - disturbance of the balance of normal flora (ATB and colon microflora) or unusual anatomic sites (perforation of colon - anaerobic peritonitis)
- Exogenous bacteria - wound contamination
- Exudate, pus, abscess
- Sampling and transport in anaerobic conditions - aspirate in the syringe with the needle capped or protected with rubber. Surgical sampling can be required. Never use the dry swab!

Anaerobic bacteria

can be divided into:

- strict anaerobes that can not grow in the presence of more than 0.5% oxygen
- moderate anaerobic bacteria that are able of growing between 2 to 8% oxygen.
- Anaerobic bacteria usually do not possess catalase
- can generate superoxide dismutase which protects them from oxygen.

- Facultative anaerobes - can grow in the presence or absence of oxygen
- Obtain energy by both respiration and fermentation
- Oxygen not toxic, some use nitrate (NO_3^-) or sulphate (SO_4^{2-}) as a terminal electron acceptor under anaerobic conditions


- Obligate (strict) anaerobes - oxygen is toxic to these organisms, do not use oxygen as terminal electron acceptor
- e.g Clostridia, Bacteriodes etc....

- **Microaerophilic organisms - require low levels of oxygen for growth, but cannot tolerate the levels present in the atmosphere**
- **Aerotolerant Anaerobes: Metabolism is anaerobic but they are unaffected by the presence of oxygen.**

Environment			
Group	Aerobic	Anaerobic	O₂ Effect
Obligate Aerobe	Growth	No growth	Required (utilized for aerobic respiration)
Microaerophile	Growth if level not too high	No growth	Required but at levels below 0.2 atm
Obligate Anaerobe	No growth	Growth	Toxic
Facultative Anaerobe (Facultative Aerobe)	Growth	Growth	Not required for growth but utilized when available
Aerotolerant Anaerobe	Growth	Growth	Not required and not utilized

Bacteria in relation to oxygen

- Obligate aerobic -require gaseous oxygen, cannot grow without oxygen
- **Obligate anaerobic - oxygen is toxic for them. Methabolic pathway used for gaining energy is fermentation with production of foul-smelling end products.**
- Facultative anaerobic - able to adapt to aerobic or anaerobic
- Microaerophilic - require reduced oxygen tension



Oxygen Toxicity

- Oxygen is used by aerobic and facultatively anaerobic organisms as its strong oxidising ability makes it an excellent electron acceptor
- During the stepwise reduction of oxygen, which takes place in respiration toxic and highly reactive intermediates are produced **reactive oxygen species (ROS)**.

The clinically important anaerobes

- 1. Gram-negative rods (*Bacteroides*, *Prevotella*, *Porphyromonas*, *Fusobacterium*,
- 2. Gram-positive cocci (primarily *Peptostreptococcus* spp.);
- 3. Gram-positive spore-forming (*Clostridium* spp.) and nonspore-forming bacilli (*Actinomyces*, *Propionibacterium*, *Eubacterium*, *Lactobacillus* and *Bifidobacterium* spp.);
- 4. Gram-negative cocci (mainly *Veillonella*)

Anaerobic environments

- Anaerobic environments (low reduction potential) include:
- Sediments of lakes, rivers and oceans; flooded soils, intestinal tract of animals; oral cavity of animals, ...
- Anaerobes also important in some infections, e.g. *C. tetanii* and *C. perfringens* important in deep puncture wound infections


- Culture of anaerobes is extremely difficult due to the need to exclude oxygen, slow growth and complex growth requirements

Culture methods

- Most common adaptation of media is the addition of a reducing agent, e.g. thioglycollate, cysteine
- Acts to reduce the oxygen to water, brings down the redox potential -300mV or less.
- Can add a redox indicator such as rezazurin, pink in the presence of oxygen - colourless in its absence


Redox potential

+500 mV



- 300 mV

OBLIGATE AEROBE MICRO AEROPHILIC FACUL-TATIVE OBLIGATE ANAEROBE



BACTERIAL GROWTH

Culture media


- Blood agar
- Blood agar - enriched with growth factors and antibiotics (selective, kanamycin, neomycin, aminoglycosides - to inhibit most aerobic and facultative anaerobic bacteria)
- VL agar, VF bouillon
- Should be stored in refrigerator in plastic bag - decrease the solution of oxygen in the agar - use freshly prepared media

Establishing of anaerobic conditions







- Candle jar - microaerophilic, facultative anaerobic
- Anaerobic jar with atmosphere generator (foil envelopes): release hydrogen and carbon dioxide -after addition of water.
- Biological method - cultivation with *Serratia marcescens*
- Tubes of broth cultures (thioglycolate, cooked meat) - their formulations contain reducing substances - anaerobic conditions at the lower part - regeneration - 15 minutes in boiling water bath

Innoculation, isolation

- Innoculation without delay
- Innoculation of anaerobic and aerobic culture (anaerobic jar, candle jar or anaerobic incubator)
- 48hrs
- Often mixed infection
- Comparison of aerobic and anaerobic culture
- Gram staining: principal is the presence of anaerobic infection G+coccus, G+rod, G-coccus, G- rod, sporulation




Algorithm


	anaerobic jar	aerobic incubator
• anaerobic		
• facultatively anaerobic		
• aerobic		



P. aeruginosa Strict aerobe




Enterococcus Facultative
Grows aerobic or anaerobic.




Bacteriodes fragilis

Anaerobic jar



(a)



(b)

Biological method

- Half of medium in the Petri's dish is inoculated with the tested sample, the other with *Serratia marcescens* - bacteria able to produce anaerobic by the oxygen. sealed with and introduced to the not oxygen - free incubator environment consumption of Petri dishes is the wax or parafin

